ENDOCRAFT

P.O. Box 8672, 25 New Chardon Street, Boston, MA 02114

Zeitels Universal Modular Glottiscope System

Sistema de Glotiscopio modular universal Zeitels

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**Zeitels Universal Modular Glottiscope System
Instructions For Use**

Limited Warranty

Endocraft warrants the original purchaser that the equipment shall be free from defects in material and workmanship for a period of one [1] year from the date of shipment. Endocraft's obligation under this warranty shall be limited to repair or replacement, at the option of Endocraft LLC. The above warranty is contingent upon normal usage and does not cover products that have been modified without Endocraft's approval or which have been subjected to unusual or unrecommended physical or chemical stress.

Notices

No part of this document may be reproduced, stored in a retrieval system, or transmitted, in any form or means -- electronic, mechanical, photocopying, recording, or otherwise -- without the expressed, written consent of Endocraft. Additional copies of this document can be ordered from Endocraft.

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MediMark® Europe, 11 Rue E. Zola, 38100 Grenoble, France

ENGLISH:**Description:**

The Zeitels Universal Modular [UM] Glottiscope system is a laryngoscope that has been designed and patented to provide improved endoscopic exposure, which enhances surgery of the true vocal folds and vestibular folds.

Indications:

The Zeitels UM Glottiscope system is indicated for endoscopic diagnosis and treatment of abnormalities of the true vocal folds and/or vestibular folds.

Assembly Instructions: (see figure on page 4)

The UM Glottiscope will be assembled in the package with the medium spatula-tube assembly [UMG101-7] attached to the universal handle and the slide assembly base plate [UMG105-SA]. The slide assembly base plate will have the detachable light cannula [UMG105-LCN] and detachable suction cannula [UMG105-SCN] inserted in position. To attach either of the other spatula-tube assemblies to the universal handle, the medium tube assembly must be detached. This is done by first separating the base-plate. The base-plate is grasped along the curved proximal sides, while holding the distal end of the spatula tube. It will separate with steady pressure. Then the spatula-tube and the vertical handle are compressed against an internal spring. This will allow a 90 degree clockwise rotation of the spatula-tube [the vertical handle is the axis of rotation]. The spatula tube can then be withdrawn and separated from the vertical handle. To attach a new spatula-tube, the process is reversed. The horizontal handle limb [UMG103-HH] is detached from the vertical handle limb [UMG102-VH] by pulling on the spring bolt and separating the two pieces. The handle extension [UMG104-HE] is detached from the horizontal handle limb by separating the two pieces with steady pressure. Attaching the components of the handle is done by reversing the aforementioned steps.

Instructions For Use:

1. Position the patient as necessary for direct laryngoscopy.
2. Place protective guards on the teeth and/or gums of the alveolar ridges.
3. The UM Glottiscope should be assembled with the desired examining spatula and without the extension limb.
4. Attach the light cable to the UM Glottiscope light carrier.
5. Place the examining tube in the oral cavity with the apex of the arch perpendicular to the teeth and advance the tube along the floor of mouth to the juncture of the oropharynx.

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ENGLISH

6. After entering the oropharynx, rotate the examining tube 90 degrees so that the apex of the arch is situated along the floor of mouth and advance the tube into the lower oropharynx.
7. Identify the tip of the epiglottis and gently advance the tip of the examining tube along the laryngeal surface of the epiglottis until the vestibular folds and/or the true vocal folds are well visualized.
8. Stabilize the UM Glottiscope in a routine fashion with the Zeitels Suspension Gallows or standard laryngoscope holder. If the suspension gallows is used, it will require that the extension limb be inserted into the horizontal handle limb.
9. Ensure that the proximal connector of the light cannula is not in contact or near the skin of the patient. A wet gauze should be placed between the light cannula and the patient's skin.
10. The light cannula can be removed completely once the microscope light is in position.
11. If a laser is used to resect tissue, the suction cannula can be placed in the base plate and connected to a vacuum source to assist with removal of the laser plume. The suction cannula is not intended for evacuation of secretions, blood, or other fluids.

For further information: Operative Techniques in Otolaryngology – Head & Neck Surgery 1999; 9:192-195

Light Guide Cannula/Fiber Optic Cable Instructions:

Intended Use:

The light guide cannula is designed to provide illumination to a surgical site by efficiently transferring light from a fiber optic light source to the laryngoscope. The light guide may be used with compatible halogen, metal halide, xenon and LED light sources.

NOTE: 300 watt or greater xenon light source should have at least 90% IR filtering to prevent damage to light guides.

NOTE: Careful matching of fiber optic cable bundle size to instrument will ensure maximum light transmission and minimize end fitting temperature.

NOTE: The light guide is provided non-sterile and must be sterilized before use. See cleaning and sterilization instructions

Instructions for Use:

1. Verify mating connections are compatible with light guide.
2. Attach light guide to the light cable and the light source and

Note: Verify all connections to light source and instruments are firm before use.

WARNING: Light emitted from this light guide contains

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significant energy. When connected to operating light source, unintended contact with the light to skin or mucosa from a light guide not attached to a device may result in burns. Do not allow distal end to directly contact tissue or other heat sensitive materials. Do not place on or near drapes, gowns or other textiles. Do not contact patient directly. Do not use for transillumination.

The temperature of the end fittings is affected by the light source and scope. These surfaces may exceed 50°C during operation. Allow light source end fitting to cool, or avoid contact after use.

CAUTION: The fiber optic cable is a precision optical instrument. Do not cut, stretch, kink, puncture or otherwise damage or alter the cable. The above will reduce light transmission and/or effect user *and patient safety*. *Discard or repair if damaged.*

ENGLISH:

How Supplied:

The Zeitels UM Glottiscope system consists of a universal handle, three examining spatula-tubes and a base-plate with detachable light and suction cannula. The universal handle is comprised of three attachable/detachable components [vertical limb, horizontal limb, extension limb]. The three examining spatula-tubes vary in size and are attachable/detachable to the vertical handle limb and the base plate.

Cleaning/Disinfecting/Sterilizing the Light Guide Cannula and Light Cable:

CLEANING:

The light guide cannula may be cleaned using a mild soap or non-oil cleaner.

CAUTION: DO NOT use synthetic detergents or oil-based soaps. The petroleum components of these soaps may be absorbed by the silicone rubber components and may leach out during use to cause a tissue reaction.

CAUTION: Avoid scratching glass fibers at ends of light guide. Damage to fibers may reduce light transmission.

NOTE: Use a mild cleaning solution with a pH range of 5 to 9.

1. Clean thoroughly using a soft-bristled brush in a lukewarm water-soap solution to remove any possible contamination.
2. Rinse thoroughly in lukewarm water.
3. Rinse thoroughly in distilled water.
4. Allow to air dry or re-sterilize, per facility policy.

STERILIZATION:

Autoclaving and Steam Sterilization:

CAUTION: DO NOT use a HIVAC sterilizer. It may cause the silicone tubing to bubble excessively.

CAUTION: Remove the fiberoptic cable from the sterilizer and allow to cool slowly to room temperature.

CAUTION: DO NOT immerse or rinse hot fiberoptic cables in cold water or liquid. This may cause glass breakage which may reduce light transmission.

Gravity Displacement Steam Sterilization:

- o Gravity displacement sterilizer, 121°C
- o Cable wrapped in an instrument tray or fully perforated sterilization box.
- o 30 minute exposure time
- o 20 minute dry time

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- o Gravity displacement sterilizer, 134°C
- o Cable unwrapped in an instrument tray or fully perforated sterilization box.
- o 4 minute exposure time
- o 20 minute dry time

Sterrad ®:

- o Use the manufacturer’s instruction for operation.
- o Do not sterilize in sterilization pouches. Doing so may result in damage to device.

Ethylene Oxide (EtO):

- o Wrapped
- o 100% EtO (725-750mg/l)
- o 1 hour 45 min exposure at 54°C
- o 12 hour aeration at 54°C

Sterilization of the Glottiscope:

Immediately after use, wipe off debris of the Glottiscope with a moistened gauze. Prior to sterilization, the suction cannula should be irrigated with at least 30 cc of saline or sterile water until it is clear. Endocraft recommends 100% ethylene oxide (EtO) gas sterilization for all components. Sterilization parameters are specified below:

EtO Sterilization (100% EtO):

Concentration:	725 mg/L (± 25 mg/L)
Temperature:	130°F (10°F/-5°F) / 54.4°C
Exposure Time:	2 Hours (+10 minutes)
Humidity:	50% (+30%/-10%)

Steam Sterilization*:

Temperature:	270 F
Exposure Time:	4 minutes
Cycle:	Pre-vacuum

*** CAUTION: THE ONLY COMPONENT OF THE UM GLOTTISCOPE THAT CANNOT BE STEAM STERILIZED IS THE LIGHT CANNULA.**

ENGLISH:

Contents

Universal Modular Glottiscope System UMG100 consists of the following:

1. **1 UMG101-9:** Large Spatula-Tube Assembly
2. **1 UMG101-7:** Medium Spatula-Tube Assembly
3. **1 UMG101-5:** Small Spatula-Tube Assembly
4. **1 UMG102-VH:** Vertical Handle Limb
5. **1 UMG103-HH:** Horizontal Handle Limb
6. **1 UMG104-HE:** Handle Extension
7. **1 UMG-LC:** Light Cable
8. **1 UMG105-SA:** Slide Assembly Base-Plate, which consists of the following:
 - a. **1 UMG105-BP:** Base Plate
 - b. **1 UMG105-LCN:** Light Cannula
 - c. **1 UMG105-SCN:** Suction Cannula

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¹ STERRAD ist eine eingetragene Marke von Johnson & Johnson Medical Inc.