# ENDOCRAFT

P.O. Box 8672, 25 New Chardon Street, Boston, MA 02114

## **Zeitels Laryngeal Suspension Gallows**

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
**ENGLISH:****Endocraft****Zeitels Universal Laryngeal Suspension Instructions  
For Use****Limited Warranty**

Endocraft warrants the original purchaser that the equipment shall be free from defects in material and workmanship for a period of one [1] year from the date of shipment. Endocraft's obligation under this warranty shall be limited to repair or replacement, at the option of Endocraft. The above warranty is contingent upon normal usage and does not cover products that have been modified without Endocraft's approval or which have been subjected to unusual or unrecommended physical or chemical stress.

**Notices**

No part of this document may be reproduced, stored in a retrieval system, or transmitted, in any form or means -- electronic, mechanical, photocopying, recording, or otherwise -- without the expressed, written consent of Endocraft. Additional copies of this document can be ordered from Endocraft.

**All inquiries, repairs, and/or returns should be directed  
to the address below:****Manufacturer****Endocraft, LLC.**

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MediMark® Europe, 11 Rue E. Zola, 38100 Grenoble, France

**ENGLISH:****Description:**

The Zeitels Laryngeal Suspension Gallows has been designed primarily to attach to the Zeitels Universal Modular Glottiscope (UMG-100A) and thereby provide elevated-vector suspension to enhance surgery of the true vocal folds and vestibular folds. It can also be used to support selected oropharyngeal mouth gags.

**Indications:**

The Zeitels Laryngeal Suspension Gallows is used to maintain the position of the UM Glottiscope for endoscopic diagnosis and treatment of abnormalities of the true vocal folds and/or vestibular folds. It can also be used to maintain the position of a number of standard oropharyngeal mouth gags.

**Unpacking and Set-up Instructions: (see exploded view of device on page 5)**

The Zeitels Laryngeal Suspension Gallows will be shipped in a collapsed position. To unpack it, ready the device for use, and attach it to the instrument rail of the operating table you must first:

1. REMOVE the clamp block from its packed position by loosening handle #2 and sliding it off the shaft.
2. ATTACH the clamp block (hang it) on the bed notch side up. Push it flat against the instrument rail and re-tighten clamping handle #1 until it locks in place.
3. PLACE the long shaft into the lumen of the clamp block it was just removed from, and advance it in until the cross pin on the shaft comes into contact with the saddle notch. Ensure that the cross pin and the top of the clamp block are in contact. Tighten handle #2.
4. RAISE the lifting mechanism into place by loosening handle #3 and swinging it up to a position 10 – 40 degrees from vertical (1-2 o'clock). After making sure the serrated teeth mesh with each other, tighten handle #3 to lock it solidly in the desired position.
5. DROPPING THE LIFTING ARM is accomplished by turning handle #4 counter clockwise till the arm is horizontal.
6. ADJUSTING THE EXTENSION SLIDE is accomplished by loosening thumbscrew #5, sliding the arm to the desired position, and then re-tightening it.
7. SETTING THE HOOK ANGLE will be accomplished by loosening & tightening black handle #6, and aligning the assembly with the laryngoscope.

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**ENGLISH:**

Handle #6 is unique since it provides the ability to re-position its angle without loosening or tightening the indwelling screw. To accomplish this, lift up the black handle and rotate it into the desired position. Upon releasing it, you may need to rotate it a small amount till it snaps back down into its new position.

**Cleaning/Disinfection:**

This device is designed to be wiped with alcohol or a comparable operating room disinfectant/equipment cleaner.

**DO NOT AUTOCLAVE THE LARYNGEAL SUSPENSION GALLOWS.**

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**ENGLISH:****Instructions For Use:**

1. Position the patient with the laryngoscope in place as necessary for direct laryngoscopy.
2. Loosening handle #1 will open the clamping jaws of the clamp block. Hang the clamp block in position on the operating room bed instrument rail in proximity to the patient's shoulder. The exact position must be individualized based on the angle of the laryngoscope and each patient's anatomy. Tighten handle #1 firmly.
3. Position the shaft of the suspension in the lumen of the clamp. Ensure that the cross pin is in contact with the top of the clamp block. Tighten handle #2 firmly.
4. Position the lifting and sliding extension arms, such that the half round cutout in the hook arm is perpendicular to the ridged horizontal extension handle of the laryngoscope. This should place the lifting mechanism approximately 10-40 degrees from the vertical. Handle #3 is tightened to maintain lifting mechanism position.
5. Then, choose a ridged slot on the extension piece of the Universal Modular Glottiscope horizontal handle into which to place the half circle of the hook arm and slowly turn handle #4 in a clockwise fashion to gently lift the laryngoscope and the patient. The suspension device will maintain its position once in place.

*For further information: Operative Techniques in Otolaryngology – Head & Neck Surgery 1999; 9:192-195*

**How Supplied:**

The Zeitels Laryngeal Suspension Gallows consists of a clamp block with two handles, a vertical shaft and lifting mechanism with handle, a sliding extension with a thumbscrew, and a hook arm with a tightening lever.

**Special Handling Instructions and \*Warnings\*:**

- ☐ **Instrumental use of any suspension device & laryngoscope can cause cervical spine trauma, dental trauma, or mucosal trauma of the oral cavity, oropharynx, hypopharynx and larynx if misused. The examiner must position the laryngoscope with care and exert lifting forces slowly, carefully monitoring the ongoing effect of the procedure on the patient to ensure that damage is not being done. Any laryngoscopy can cause temporary numbness of the tongue from steady pressure while the laryngoscope is positioned with a holder or suspension gallows.**

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**ENGLISH:**

- ❑ There are patients who cannot be exposed adequately by any direct laryngoscope and therefore undue force should not be exerted.
- ❑ **Extra care** should be observed when extending the head and neck of any patient with a rigid or abnormal cervical spine, especially if there is known cervical spine disease. Instrumental use of any laryngeal suspension device can cause trauma to the cervical spine if excessive force is used and especially in those with a variety of spine problems such as degeneration or prior trauma. The use of a direct laryngoscope with a suspension device can cause trauma to the temporal-mandibular joint, teeth or mucosal trauma to the oral cavity, oropharynx, hypopharynx or larynx if misused.
- ❑ Interchange of any components with another company's instrumentation such as attaching an alternative laryngoscope to the Zeitels Laryngeal Suspension Gallows is generally not suggested and must be assessed for safety and efficacy by the treating physician.

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**ENGLISH:**Contents

Zeitels Laryngeal Suspension Gallows ZVG-100A consists of the following:

1. ZVG-101-HA Hook arm assembly
2. ZVG-102-SA Slide arm assembly
3. ZVG-104-GB Gear box assembly
4. ZVG 104-LA Lifting Arm assembly
5. ZVG-105-SA Swing arm assembly
6. ZVG-106-VS Vertical Shaft
7. ZVG-107-CB: Clamp Block Assembly

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